


# Preprocessing fMRI data using BIDS and fMRIPrep

Neal Morton

# Preprocessing fMRI data



- Preprocessing tools
- BIDS format
- fMRIPrep

# Preprocessing tools

# Preprocessing tools

## Limitations

- Many individual tools: AFNI, SPM, FSL, FreeSurfer, ANTS...
- Installation is complicated and idiosyncratic
- The interface and language used for each package varies
- For optimal preprocessing, need to combine tools that were not designed to work together (e.g., FSL and ANTS)

# Preprocessing tools


## Improvements

- Interoperability made possible by the NIfTI format (toolbox-specific formats remain in FreeSurfer and AFNI, however)
- To further help bridge tools, Nipype was developed
- Docker containers have further helped ease use of multiple software packages together





# Nipype

- Describes "workflows" as a graph of inputs, operations, and outputs
- A replacement for shell scripts
- Distributes processes over available cores and manages memory usage
- Workflows are available and under active development for multiple modalities
  - Spatial preprocessing
  - Susceptibility distortion correction
  - fMRIPrep



# Docker containers


- Containers allow apps, including their dependencies (e.g., FSL, ANTs), to be easily installed and run
- Containers are created from a Dockerfile that specifies how to install everything needed to run the app
- Containers can be downloaded from Docker Hub
- Docker requires root permissions, which doesn't work on a cluster; enter Singularity
- Singularity works on TACC resources, and a number of containers are installed to run software used in the lab





# fMRIPrep

- Uses Nipype to implement "minimal preprocessing" of fMRI datasets
- Includes anatomical segmentation, slice-timing correction, motion correction, distortion correction, and registration to template space
- Idea is to run a common set of preprocessing for any analysis (univariate, multivariate, connectivity, etc.)
- Uses many different tools under the hood, which are all packaged in a Docker container
- Requires data to be in BIDS format



# Preprocessing tools


- See the lab wiki for links to documentation for each of these tools
  - <https://github.com/prestonlab/wiki/wiki/Running-basic-fMRI-preprocessing-using-BIDS-and-fMRIPrep>

# BIDS format

# BIDS

## Brain imaging data structure

- Can streamline software development and improve usability with a standard data structure
- Put your data into BIDS format once, then easily run any "BIDS app"
- Data are stored in standard formats (NIfTI, TSV), while metadata are stored in "sidecar" JSON files
- Encourages complete documentation, including scanning parameters and behavioral task events
- Most metadata are extracted automatically





# Converting task data

Example:

```
└─ sub-control01/
   └─ func/
      ├── sub-control01_task-emoface_events.tsv
      └── sub-control01_task-emoface_events.json
```

Example of the content of the TSV file:

```
onset duration trial_type identifier database response_time
1.2 0.6 afraid AF01AFAF kdef 1.435
5.6 0.6 angry AM01AFAN kdef 1.739
5.6 0.6 sad AF01ANSA kdef 1.739
```

The `trial_type` and `identifier` columns from the `events.tsv` files might be described in the accompanying JSON sidecar as follows:

```
{
  "trial_type": {
    "LongName": "Emotion image type",
    "Description": "Type of emotional face from Karolinska database that is displayed",
    "Levels": {
      "afraid": "A face showing fear is displayed",
      "angry": "A face showing anger is displayed",
      "sad": "A face showing sadness is displayed"
    }
  },
  "identifier": {
    "LongName": "Karolinska (KDEF) database identifier",
    "Description": "ID from KDEF database used to identify the displayed image"
  }
}
```

Note that all other columns SHOULD also be described but are omitted for the sake of brevity.

- The BIDS specification includes a simple, flexible format for describing experiment events
- TSV (tab-separated value) file for each run. Can write using any language (Python, R, Matlab, all have good tools for working with tabular data)
- Generally should define onset, duration, and trial\_type fields. Add any other fields you need for neural or behavioral analysis
- Fields are described in a JSON sidecar file (can just have one per task, in the main directory of the BIDS structure)
- Can write directly to events format when collecting data, or write a script to create events files after the fact (see Tesser and Bender for examples)



# HeuDiConv

## Heuristic Dicom Conversion

- HeuDiConv can be configured to convert basically any fMRI dataset to BIDS format
- First, extract information from the DICOM files to get features of each scan
- For each dataset, need to write a Python "heuristic file" that uses those features to put each scan in the correct place
- Once the heuristic is determined, can use HeuDiConv to run the dcm2niix program to convert to NIfTI format
- The main HeuDiConv page doesn't really explain how to use it in practical terms
  - See the Stanford tutorial for a walkthrough

dicominfo\_ses-001 — Saved to my Mac

View View

Wrap Text Merge & Center

General

\$ %

0.00 0.00

G	H	I	J	K	L	M
dim1	dim2	dim3	dim4	TR	TE	protocol_name
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra

```
import os
def create_key(template, outtype=('nii.gz',), annotation_classes=None):
    if template is None or not template:
        raise ValueError('Template must be a valid format string')
    return template, outtype, annotation_classes
def infotodict(seqinfo):
    """Heuristic evaluator for determining which runs belong where
    allowed template fields - follow python string module:
    item: index within category
    subject: participant id
    seqitem: run number during scanning
    subindex: sub index within group
    """
    t1w = create_key('sub-{subject}/{session}/anat/sub-{subject}_{session}_run-00{item:01d}_T1w')
    info = {t1w: []}
    for idx, s in enumerate(seqinfo):
        if (s.dim1 == 320) and (s.dim2 == 320) and ('t1_fl2d_tra' in s.protocol_name):
            info[t1w].append(s.series_id)
    return info
```



# HeuDiConv

## Recommended workflow

- Locally, install just the Python heudiconv package (don't need dcm2niix)
- Follow the walkthrough to get a scan spreadsheet (see results in hidden `.heudiconv` directory)
- Copy the template heuristic file and customize it for your scans
- Run heudiconv with your heuristic file with no conversion (`-c none`) to see how scans will be named
- Can use the Python debugger to better understand how the code runs (add `breakpoint()` in the heuristic file)
- On the cluster, use the heuristic file to run the actual conversion using Singularity (see Tesser and Bender projects for examples)

dicominfo\_ses-001 — Saved to my Mac

View View

Wrap Text Merge & Center

General

\$ %

0.00 0.00

G	H	I	J	K	L	M
dim1	dim2	dim3	dim4	TR	TE	protocol_name
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra
320	320	10	1	0.15	3.74	t1_fl2d_tra

```
import os
def create_key(template, outtype=('nii.gz',), annotation_classes=None):
    if template is None or not template:
        raise ValueError('Template must be a valid format string')
    return template, outtype, annotation_classes
def infotodict(seqinfo):
    """Heuristic evaluator for determining which runs belong where
    allowed template fields - follow python string module:
    item: index within category
    subject: participant id
    seqitem: run number during scanning
    subindex: sub index within group
    """
    t1w = create_key('sub-{subject}/{session}/anat/sub-{subject}_{session}_run-00{item:01d}_T1w')
    info = {t1w: []}
    for idx, s in enumerate(seqinfo):
        if (s.dim1 == 320) and (s.dim2 == 320) and ('t1_fl2d_tra' in s.protocol_name):
            info[t1w].append(s.series_id)
    return info
```

# Running HeuDiConv on TACC

- On TACC clusters like Lonestar 6, can only run Singularity on compute nodes, not login nodes
- Load the tacc-singularity module
- Specify the path to a container or "image", followed by the inputs to the program
- Containers have entry points that work the same as scripts (here, same as running the heudiconv command)

```
image=$STOCKYARD2/software/images/heudiconv-0.9.0.sif
module load tacc-singularity
singularity run "$image" \
    -s "$subject" \
    -f "$heuristic" \
    -b \
    -o "$bids_dir" \
    --minmeta \
    --files "$subj_raw_dir"/*/*.dcm
```

# HeuDiConv tips

- Convert both the BOLD timeseries images and the "SBRef" (single-band reference scans)
  - fMRIPrep can use SBRef scans for improved registration between functional and anatomical scans
- Exclude any scans you won't analyze (e.g., aborted scans) at this stage. Can check the number of volumes and exclude scans with less than the expected number
- Can program exceptions in the heuristic file (e.g., if `patient_id` is 102, exclude the first coronal scan)
- Can run into problems with multi-day studies, because scan order is determined based on a number assigned within the scanning session. Must first edit this number to avoid problems (have a Python script to do this in the Bender GitHub project)



# After conversion

- After conversion, need to set the IntendedFor field for the JSON files associated with each fieldmap scan
- fMRIPrep uses this information when correcting for distortion of functional scans
- The PyBIDS package can help work with BIDS datasets. The Tesser and Bender GitHub projects use PyBIDS to add IntendedFor information to the fieldmap scans

```
"IntendedFor": [  
    "func/sub-30_task-localizer_run-1_bold.nii.gz",  
    "func/sub-30_task-localizer_run-2_bold.nii.gz",  
    "func/sub-30_task-localizer_run-3_bold.nii.gz",  
    "func/sub-30_task-localizer_run-4_bold.nii.gz",  
    "func/sub-30_task-exposure_run-1_bold.nii.gz"  
]
```




# Validating a dataset

- Can use the BIDS-Validator program to check a converted dataset
  - `singularity run /work/03206/mortonne/software/images/validator-1.7.3.sif path_to_bids_directory`
- Will flag missing or mismatching files (e.g., if participants have different scans)
- May need to adjust your heuristic file and run HeuDiConv again


fMRIprep

# fMRIprep workflow






# Preprocessing steps




# Preprocessing steps





# Organizing data

- Have three main directories:
  - sourcedata - raw DICOM files, raw behavioral files, etc.
  - rawdata - BIDS format neural and behavioral data
  - derivatives - results from analyzing BIDS data



# Running fMRIPrep

- The fMRIPrep webpage has an example script to run it using Singularity
- Takes in BIDS format scan data
- Writes intermediate files to a "work" directory (not to be confused with the TACC WORK filesystem) when running
  - \$SCRATCH a good filesystem to use for this
- In the derivatives directory, only outputs the data you request, without intermediate files
- Templates are downloaded to a TemplateFlow directory

```
# Set singularity command and mount points
SINGULARITY_CMD="singularity run \
  --cleanenv \
  -B ${BIDS_DIR}:/data \
  -B ${TEMPLATEFLOW_HOST_HOME}:${SINGULARITYENV_TEMPLATEFLOW_HOME} \
  -B ${WORKDIR}:/workdir \
  ${STOCKYARD2}/software/images/fmriprep-20.2.1.simg"

# Compose the command line
cmd="${SINGULARITY_CMD} \
  /data/rawdata \
  /data/${DERIVS_DIR} \
  participant \
  --participant-label ${subject} \
  -w /workdir/ \
  --fd-spike-threshold 0.5 \
  --dvars-spike-threshold 1.5 \
  --output-spaces MNI152NLin2009cAsym:res-native fsaverage:den-164k anat:res-native \
  -v \
  --omp-nthreads 12 \
  --nthreads 18 \
  --mem_mb 60000"
```

# Output spaces

- Preprocessed data can be written in multiple output spaces
  - Native anatomical space
  - MNI template space
    - The MNI152NLin2009cAsym template is a better version of the FSL template
  - Surface space
    - fsaverage
    - fsLR (seems to be a newer version)
- Can choose the resolution (e.g., template space but in native resolution at 1.7 mm iso)

```
# Set singularity command and mount points
SINGULARITY_CMD="singularity run \
  --cleanenv \
  -B ${BIDS_DIR}:/data \
  -B ${TEMPLATEFLOW_HOST_HOME}:${SINGULARITYENV_TEMPLATEFLOW_HOME} \
  -B ${WORKDIR}:/workdir \
  ${STOCKYARD2}/software/images/fmriprep-20.2.1.simg"

# Compose the command line
cmd="${SINGULARITY_CMD} \
  /data/rawdata \
  /data/${DERIVS_DIR} \
  participant \
  --participant-label ${subject} \
  -w /workdir/ \
  --fd-spike-threshold 0.5 \
  --dvars-spike-threshold 1.5 \
  --output-spaces MNI152NLin2009cAsym:res-native fsaverage:den-164k anat:res-native \
  -v \
  --omp-nthreads 12 \
  --nthreads 18 \
  --mem_mb 60000"
```

# fMRIPrep output

## Quality assurance

- Derivatives like fMRIPrep are formatted similarly to BIDS raw data
- Get an HTML report for quality assurance for each participant
- Download the html files and figures directories, then view the html files in a browser

```
fmriprep
├── dataset_description.json
├── desc-aparcaseg_dseg.tsv
├── desc-aseg_dseg.tsv
├── logs
│   ├── CITATION.bib
│   ├── CITATION.html
│   ├── CITATION.md
│   └── CITATION.tex
├── sub-01
│   ├── anat
│   ├── figures
│   ├── func
│   └── log
├── sub-01.html
├── sub-20
│   ├── anat
│   ├── figures
│   ├── func
│   └── log
├── sub-20.html
├── sub-25
│   ├── anat
│   ├── figures
│   ├── func
│   └── log
├── sub-25.html
├── sub-30
│   ├── anat
│   ├── figures
│   ├── func
│   └── log
└── sub-30.html
```

# fMRIprep output

## Anatomical

- FreeSurfer ROIs are in the dseg files
  - See desc-aseg\_dseg.tsv in the main directory for labels
- Spatial transformations specified in .txt and .h5 files
- Surface geometry is in the .gii files
- If no space is indicated, it is in native space
- In each space:
  - T1 scan
  - Brain mask
  - Basic segmentation (gray matter, white matter, CSF)

fmriprep/sub-01/anat

```
— sub-01_desc-aparcaseg_dseg.nii.gz
— sub-01_desc-aseg_dseg.nii.gz
— sub-01_desc-brain_mask.json
— sub-01_desc-brain_mask.nii.gz
— sub-01_desc-preproc_T1w.json
— sub-01_desc-preproc_T1w.nii.gz
— sub-01_dseg.nii.gz
— sub-01_from-fsnative_to-T1w_mode-image_xfm.txt
— sub-01_from-MNI152NLin2009cAsym_to-T1w_mode-image_xfm.h5
— sub-01_from-T1w_to-fsnative_mode-image_xfm.txt
— sub-01_from-T1w_to-MNI152NLin2009cAsym_mode-image_xfm.h5
— sub-01_hemi-L_inflated.surf.gii
— sub-01_hemi-L_midthickness.surf.gii
— sub-01_hemi-L_pial.surf.gii
— sub-01_hemi-L_smoothwm.surf.gii
— sub-01_hemi-R_inflated.surf.gii
— sub-01_hemi-R_midthickness.surf.gii
— sub-01_hemi-R_pial.surf.gii
— sub-01_hemi-R_smoothwm.surf.gii
— sub-01_label-CSF_probseg.nii.gz
— sub-01_label-GM_probseg.nii.gz
— sub-01_label-WM_probseg.nii.gz
— sub-01_run-1_from-orig_to-T1w_mode-image_xfm.txt
— sub-01_run-2_from-orig_to-T1w_mode-image_xfm.txt
— sub-01_space-MNI152NLin2009cAsym_desc-brain_mask.json
— sub-01_space-MNI152NLin2009cAsym_desc-brain_mask.nii.gz
— sub-01_space-MNI152NLin2009cAsym_desc-preproc_T1w.json
— sub-01_space-MNI152NLin2009cAsym_desc-preproc_T1w.nii.gz
— sub-01_space-MNI152NLin2009cAsym_dseg.nii.gz
— sub-01_space-MNI152NLin2009cAsym_label-CSF_probseg.nii.gz
— sub-01_space-MNI152NLin2009cAsym_label-GM_probseg.nii.gz
— sub-01_space-MNI152NLin2009cAsym_label-WM_probseg.nii.gz
```



# fMRIPrep output

## Functional

- Preprocessed BOLD data are in the desc-preproc\_bold files
- Here, written in three spaces:
  - Native anatomical (useful for ROI analysis)
  - MNI template space (useful for univariate analysis and searchlights)
  - Surface space (can be used for both individual and group analysis)
- Also have confound information (motion parameters, CSF signal etc.)

fmriprep/sub-01/func

```
— sub-01_task-ACtest_run-1_desc-confounds_timeseries.json
— sub-01_task-ACtest_run-1_desc-confounds_timeseries.tsv
— sub-01_task-ACtest_run-1_from-scanner_to-T1w_mode-image_xfm.txt
— sub-01_task-ACtest_run-1_from-T1w_to-scanner_mode-image_xfm.txt
— sub-01_task-ACtest_run-1_space-fsaverage_hemi-L_bold.func.gii
— sub-01_task-ACtest_run-1_space-fsaverage_hemi-L_bold.json
— sub-01_task-ACtest_run-1_space-fsaverage_hemi-R_bold.func.gii
— sub-01_task-ACtest_run-1_space-fsaverage_hemi-R_bold.json
— sub-01_task-ACtest_run-1_space-MNI152NLin2009cAsym_boldref.nii.gz
— sub-01_task-ACtest_run-1_space-MNI152NLin2009cAsym_desc-aparcaseg_dseg.nii.gz
— sub-01_task-ACtest_run-1_space-MNI152NLin2009cAsym_desc-aseg_dseg.nii.gz
— sub-01_task-ACtest_run-1_space-MNI152NLin2009cAsym_desc-brain_mask.json
— sub-01_task-ACtest_run-1_space-MNI152NLin2009cAsym_desc-brain_mask.nii.gz
— sub-01_task-ACtest_run-1_space-MNI152NLin2009cAsym_desc-preproc_bold.json
— sub-01_task-ACtest_run-1_space-MNI152NLin2009cAsym_desc-preproc_bold.nii.gz
— sub-01_task-ACtest_run-1_space-T1w_boldref.nii.gz
— sub-01_task-ACtest_run-1_space-T1w_desc-aparcaseg_dseg.nii.gz
— sub-01_task-ACtest_run-1_space-T1w_desc-aseg_dseg.nii.gz
— sub-01_task-ACtest_run-1_space-T1w_desc-brain_mask.json
— sub-01_task-ACtest_run-1_space-T1w_desc-brain_mask.nii.gz
— sub-01_task-ACtest_run-1_space-T1w_desc-preproc_bold.json
— sub-01_task-ACtest_run-1_space-T1w_desc-preproc_bold.nii.gz
```

# Conclusions

- BIDS format makes neuroimaging data better documented and easier to work with
- BIDS apps can run relatively simply without the user having to specify many options; Docker/Singularity containers simplify installation
- fMRIPrep is supported by an active community of developers who regularly fix bugs and add new features
- fMRIPrep output is ready for a number of analysis approaches